# APS Data Processing

#### SRJ Likith

# 1 Steps

#### 1 Read relevant files

First, the files containing the data and information that we need are read into pandas DataFrames using the read\_excel and read\_csv commands as appropriate.

```
import pandas as pd
legend_df=pd.read_excel('SAXSframevstest_final.xlsx')
meta_df=pd.read_csv('aclarke_jul21_exp_tracking.csv',low_memory=False)
```

### 2 Cleaning the data

The first cleaning step, and these are aimed mainly at the metadata file, is to strip column names of leading and trailing whitespace characters. Next, any columns that consist of all null, None, or NaN values are removed. Lastly, the Date column is converted from strings to datetime objects. This is all achieved by the clean\_df function and the helper function dt\_formatter.

```
from datetime import datetime as dt

def dt_formatter(dt_str):
    return dt.strptime(dt_str, '%a %b %d %H:%M:%S %Y)

def clean_df(in_df):
    in_df.columns=in_df.columns.str.strip()
    for c in in_df.columns:
        if in_df[c].isnull().all():
            del in_df[c]
    in_df['Date']=in_df['Date'].apply(dt_formatter)
    return in_df

meta_df=clean_df(meta_df)
```

### 3 Adding calculated column(s)

Based on the saxs fnum column, the wax\_fnum column is added, which is saxs fnum/10 rounded to the nearest integer.

```
def wax_fnum(row):
    return round(int(row['saxs fnum'])/10)

meta_df['wax fnum']=meta_df.apply(wax_fnum,axis=1)
```

#### 4 Processing

In order to achieve all the processing steps, we loop through (for loop) the legend\_df DataFrame. The idx variable contains the index/row number while the row variable takes the values of the rows themselves.

```
for idx, row in lengend_df.iterrows():
```

Note that the most following lines are all contained within this for loop and so, are indented one level (4 spaces/1 tab).

For each row, the first thing to do is to read the SAXS Frame start and SAXS Frame End, and then apply the buffer value, i.e., the limits now become SAXS Frame start-buffer\_frames and SAXS Frame end+buffer\_frames. Note that the buffer\_frames=100 line needs to be run just once, and so, should be placed before the start of the for loop.

```
buffer_frames=100 #this line needs to be placed before the for loop starts
fnum_lims=[int(i) for i in [row['SAXS Frame Start'],row['SAXS Frame End']]]
fnum_lims=[el+((-1)**(idex+1)*buffer_frames) for idex,el in enumerate(fnum_lims)]
```

Next, the sample name is read from the SampleName column and stripped of leading and trailing whitespaces. The fname variable is meant for naming the files and directories, and since this shouldn't have any / characters, these are replaced by \_ characters.

```
sample=row['SampleName'].strip()
fname=sample.replace('/','_')
```

Then, a small subset/subsection of the meta\_df is chosen that consists of the relevant columns, and only those rows with saxs fnum that lie between the limits mentioned earlier. Also, the MTS load (mm) is renamed to Force (N). Again, like the previous block of code, these lines are indented one level since they are placed within the for loop.

Next, for each row in legend\_df, a folder is created (if it doesn't already exist), and the curr\_meta DataFrame is written to that directory. Note that the file\_ops function is defined before the loop. The variable p helps point to the directory where we want all the sections and their directories deposited. This folder is created if it doesn't already exist.

```
from pathlib import Path
p=Path('./sections_dir_final')
p.mkdir(parents=True, exist_ok=True)
def file_ops(sample_name,dest_dir=p)
target=dest_dir/sample_name
if not target.exists():
target.mkdir(parents=True,exist_ok=True)

##the following lines are contained within the for loop, and hence, indented
file_ops(fname)
curr_meta.to_csv(p/fname/(fname+'_metadata.csv'),sep=',', index=False)
```

Finally, a small readme file is deposited in each of the sub-directories, containing the following information about buffer\_frames

# 2 Appendix

```
#!/usr/bin/python3
3 from pathlib import Path
4 import pandas as pd
5 from glob import glob
6 from datetime import datetime as dt
s def dt_formatter(dt_str): #function to format the 'Date' column in
     meta_df ie dataframe containing the metadata file
     return dt.strptime(dt_str, '%a %b %d %H:%M:%S %Y')
def clean_df(in_df): #cleaning the meta_df. strips column names of
     leading and trailing whitespaces. deletes any columns where all the
     values are null, NaN, etc.
     in_df.columns=in_df.columns.str.strip()
12
     for c in in_df.columns:
          if in_df[c].isnull().all():
14
              del in_df[c]
      in_df['Date'] = in_df['Date'].apply(dt_formatter)
16
     return in_df
17
def wax_fnum(row):
     return round(int(row['saxs fnum'])/10)
p=Path('./sections_dir_final')
p.mkdir(parents=True, exist_ok=True)
def file_ops(sample_name,dest_dir=p):
      target=dest_dir/sample_name
      if not target.exists():
27
          target.mkdir(parents=True, exist_ok=True)
pd.set_option('display.max_columns', None)
1 legend_df=pd.read_excel('SAXSframesvstest_final.xlsx') #contains SAXS
     frame start, end numbers and SampleNames
32 buffer_frames=100 #number of frames added to the beginning and end of
     each section to act as a buffer for context i.e., each section,
     instead of consisting of SAXS Frame Start to SAXS Frame End, will
     consist of (SAXS Frame Start - buffer) to (SAXS Frame End + buffer)
33 readme_str='buffer_frames=100 #number of frames added to the beginning
     and end of each section to act as a buffer for context i.e., each
     section, instead of consisting of SAXS Frame Start to SAXS Frame End,
     will consist of (SAXS Frame Start - buffer) to (SAXS Frame End +
     buffer)'
meta_df=pd.read_csv('aclarke_jul21_exp_tracking.csv',low_memory=False)
meta_df = clean_df (meta_df)
meta_df['wax fnum']=meta_df.apply(wax_fnum,axis=1)
for idx, row in legend_df.iterrows():
```

```
fnum_lims=[int(i) for i in [row['SAXS Frame Start'],row['SAXS Frame
     End ']]]
      fnum_lims=[el+((-1)**(idex+1)*buffer_frames) for idex,el in
41
     enumerate(fnum_lims)]#subtracting buffer_frames from frame start and
     adding it to frame end
      sample=row['SampleName'].strip()
      fname=sample.replace('/','_')
43
      print(sample)
      curr_meta=meta_df[['Date', 'MTS load (mm)', 'Furnace T1 (C)',
45
                          'saxs fnum', 'wax fnum', 'MTS crosshead
     (mm)']][(meta_df['saxs fnum']>=fnum_lims[0]) & (meta_df['saxs
     fnum'] <= fnum_lims[1])]</pre>
      curr_meta.rename({'MTS load (mm)':'Force (N)'}, axis=1, inplace=True)
47
      file_ops(fname)
48
      curr_meta.to_csv(p/fname/(fname+'_metadata.csv'),sep=',',index=False)
      with open(p/fname/'readme.txt','w') as f:
50
          f.writelines(readme_str)
```

Listing 1: processor\_final.py